

BEEA[®]SKIN

BEEA[®]SKIN

The future of skincare, in your hands.

BEEA SKIN

Brand Philosophy

Beea Skin brings professional skincare into the home through advanced, science-driven beauty devices focused on anti-aging. Using **RF**, **EMS**, and **LED technologies**, the brand stimulates collagen and revitalizes skin, while offering a holistic approach for face, body, and hair. With upcoming complementary skincare products, Beea Skin promotes confidence, radiance, and everyday wellness.

Beea Skin seamlessly integrates beauty and self-care into daily routines, guided by innovation, effectiveness, and user empowerment—transforming skincare into an inspiring, accessible experience.

Brand Promise

We empower you with state-of-the-art skincare devices that deliver professional results at home. Driven by innovation and science, our solutions are tailored to your unique needs. As your trusted partner, we are expanding into complementary skincare products to create a complete, results-focused beauty experience.





Beea Skin Longevity

The Strategic Standard for 2026

In 2026, BEEA SKIN sets a higher bar: **skin longevity**. Not a promise of a quick effect — but a system built for **long-term skin stability**. Powered by **BEEA SKIN devices** and a structured at-home methodology, skin longevity becomes a repeatable standard — not a one-time result.

Because the real goal is not “better today.”
The real goal is skin that, over time, becomes:

- more stable
- more even
- more resilient
- more consistent in quality and appearance

Skin longevity in BEEA SKIN is a measurable methodology. Built on disciplined application, clearly defined control points, and verified progress — not assumptions.

Skin longevity is not a trend.
It's the future of skincare — executed as a standard.

Categories




CLEANSING



ANTI - AGING



HAIR CARE



BODY CARE


BEEA SKIN



AURA, BEEA SKIN, LIFTING AND ANTI-AGING DEVICE

BEEA AURA is a multifunctional facial skincare device that combines the most advanced technologies: microcurrent, EMS, electroporation, and LED light therapy. Designed for users seeking precise contouring, elasticity restoration, and pore treatment – all in one sophisticated device.

- stimulates collagen production
- reduces puffiness
- reduces redness and inflammation
- skin regeneration
- healthy and even glow
- reduces wrinkles and fine lines
- tightens and tones the skin
- clears acne



BS-2370
299,90 €



DEEP DERMAL HYDRATION BOOST

Enhances transdermal delivery of active ingredients



COLLAGEN & ELASTICITY REACTIVATION

Microcurrent + EMS stimulation activates fibroblast activity and facial muscle fibers



FACIAL MUSCLE TONING & LIFT

EMS technology stimulates underlying facial musculature & improves muscle tone and contour definition.



ACCELERATED SKIN RENEWAL

Supports cellular turnover and recovery processes



MULTI-LAYER REJUVENATION SYSTEM

Works across multiple biological levels



EXPECTED RESULTS WITH REGULAR USE*

- Increased skin hydration: +30–40%
- Improved elasticity: +18–25%
- Enhanced facial muscle tone: visible contour refinement
- Reduction in appearance of fine lines within 4–8 weeks
- Improved overall radiance and skin density

*Results based on internal performance evaluations and typical response to combined microcurrent/EMS electroporation technologies.


BEEA[®]SKIN



SCULPTA, BEEA SKIN, FACE LIFT DEVICE

The device with a 3D curved head adapts to the contours of the face – especially the chin and neck – to firm and tighten the skin. It combines LED light therapy (red, green, blue) with EMS technology for improved skin texture. Three functions in one device: lifting, skin tone evening, and soothing – your complete beauty routine in just 6 minutes a day.

- boosts collagen
- reduces puffiness
- calms irritation
- clears acne
- healthy, even glow
- smooths lines and wrinkles
- tightens and tones
- skin regeneration




EMSTECHNOLOGY



LED LIGHT



HOT THERAPY



VIBRATION TECHNOLOGY

BS-2368
129,90 €

6-MINUTE DAILY LIFT PROTOCOL

Muscle stimulation

LED therapy

Thermal activation

Vibration-assisted circulation

Efficient. Targeted. Structured.

EMS MUSCLE ACTIVATION & TONING

Supports visibly lifted appearance with consistent use



COLLAGEN & TEXTURE SUPPORT

Red and green LED light therapy promotes collagen activity



THERMAL STIMULATION FOR FIRMNESS

*Controlled heat therapy enhances microcirculation
and supports collagen remodeling*



LOWER FACE CONTOUR DEFINITION

*Ergonomic 3D curved head precisely
adapts to the jawline and neck*




BEEA[®] SKIN



REMS, BEEA SKIN, ANTI-AGING DEVICE

The REMS beauty device combines three advanced technologies for skin rejuvenation and care. It is easy to use, suitable for all skin types, and delivers visible results quickly. Achieve firmer, more radiant skin from the comfort of your home.

- rejuvenated, tightened skin
- improved skin tone
- lifting effect
- glow effect
- firm and radiant skin
- less inflammation and irritation
- evened-out complexion
- better circulation



EMSTECHNOLOGY



LED LIGHT



RF TECHNOLOGY



INFRARED LIGHT
THERAPY

BS-2382
199,90 €



2 MHz RADIOFREQUENCY SKIN TIGHTENING

Delivers controlled dermal heating to stimulate collagen

RED LIGHT COLLAGEN ACTIVATION

Supports fibroblast activity and collagen synthesis



TRIPLE ENERGY REJUVENATION SYSTEM

*combination of RF, Red LED, and Infrared
Light works across multiple skin layers*



LIFT. FIRM. REVITALIZE.

Clinically inspired energy-based technology



INFRARED DEEP TISSUE STIMULATION

*support circulation, accelerate
recovery processes,*

**Helena, 52 – results after 1 month
- REMS device**

- excess skin and under-eye puffiness reduced
- improved skin tone across the entire face
- visibly noticeable anti-aging effect
- improved skin texture
- more even complexion, purified and revitalized appearance
- restored tone to skin and facial muscles
- less pronounced wrinkles and fine lines
- lips appear fuller with better definition
- skin appears firmer, more elastic, and hydrated



before

after



before

after

BEEA SKIN



FACE GYM, BEEA SKIN, ANTI-AGING DEVICE FOR FACE, EYES & LIPS

This versatile device combines three advanced technologies in one: microcurrent (MC), EMS (electrical muscle stimulation), and LED therapy. With daily use, it trains facial muscles, improves their tone, and defines contours. Featuring an elegant design and suitable for everyday use, this beauty device will transform and elevate your skincare routine. Face Gym 3-in-1 comes with 3 interchangeable heads, each designed for targeted treatments.

reduces acne and inflammation

reduces puffiness and dark circles

evens out skin tone


boosts collagen production

lifts and sculpts the face


improves circulation and skin tone

smooths wrinkles and fine lines

speeds up skin regeneration



MICROCURENT
TECHNOLOGY



EMS TECHNOLOGY



LED LIGHT

TB-1767
199,90 €

Martina, 46, results after 2 months
- FaceGym device, REMS

- skin appears firmer, more elastic, and visibly revitalized
- pores appear tighter and the skin looks more vibrant
- cheekbones appear lifted and fuller
- under-eye area looks brighter and better circulated
- skin looks fresher and more even
- nasolabial folds are softened
- fine lines at the outer corners of the lips are reduced
- barcode lines (smoker's lines) appear less pronounced
- under-eye area looks brighter and better circulated



before

after



before

after

BEEA SKIN



GLAM, BEEA SKIN, LED FACE MASK

The BEEA Skin LED mask provides safe and effective at-home skincare, improves tone and texture, and enables professional rejuvenation without UV radiation.

- | | |
|----------------------------|----------------------|
| skin regeneration | skin firming |
| cleaner and smoother skin | fewer signs of aging |
| soothes and hydrates | evens out complexion |
| less redness & sensitivity | less clogged pores |



LED LIGHT



RED LIGHT



BLUE LIGHT



YELLOW LIGHT

BS-2386
199,90 €

ANTI-INFLAMMATORY & BARRIER SUPPORT

Yellow light at 590 nm helps calm erythema, reduce sensitivity, and enhance skin recovery



HIGH-INTENSITY PROFESSIONAL IRRADIANCE

Red: 100 mW/cm² Yellow: 70 mW/cm² Blue: 550 mW/cm²



COLLAGEN & FIRMNESS ACTIVATION

Red light at 625 nm penetrates dermal layers to stimulate fibroblast activity



TARGETED ACNE CONTROL TECHNOLOGY

Blue light at 415 nm with high irradiance (550 mW/cm²) supports antibacterial action against acne-related bacteria while helping regulate breakouts



PRECISION WAVELENGTH THERAPY

*Red Light – 625 ± 10 nm → collagen stimulation & wrinkle reduction
Yellow Light – 590 ± 10 nm → redness calming & microcirculation support
Blue Light – 415 ± 10 nm → acne-targeting antibacterial action*



BEEA SKIN




VITA, BEEA SKIN, ANTI-AGING HAIR COMB

Integrated LLLT technology (Low Level Laser Therapy), combined with other technologies, makes this device ideal for anyone struggling with hair loss and looking to prevent it. LLLT controls sebum secretion on the scalp, prevents hair loss, and strengthens and activates existing hair follicles.

- | | |
|-----------------------------------|----------------------------------|
| smooth, shiny hair without static | absorption of active ingredients |
| detangling without pulling | less oiliness and dandruff |
| improved circulation | stimulates hair growth |
| stronger hair roots | scalp revitalization |
| soothes irritation | healthy-looking hair |



LED LIGHT



RED LIGHT



BLUE LIGHT



YELLOW LIGHT

BS-2343
199,90 €



**STRONGER ROOTS.
HEALTHIER GROWTH.**

Stimulates the scalp




**CLINICALLY-
INSPIRED HAIR
ACTIVATION**

Advanced LLLT technology




**MULTI-TECH SCALP
REJUVENATION
SYSTEM**

LLLT, EMS, LED light



**ENHANCED ABSORPTION
& SHINE BOOST**

Improves penetration of active ingredients



**SCALP OIL &
DANDRUFF BALANCE**

Helps regulate sebum production,



before



after

VITA brush

- scalp revitalization with red LED light
- oil regulation with blue LED light
- improved microcirculation (sonic vibration)
- extended hair growth phase (laser stimulation)
- negative ions for smoother hair, reduced static, and enhanced shine

BEEA[®]SKIN


NEW LAUNCHES

LUMA, BEEA SKIN, LED PHOTON MASK

LUMA delivers photobiomodulation on two depths at once. 84 quad-chip LEDs cover the full face with clinically mapped wavelengths — 630 nm red for collagen, 590 nm yellow for redness, 415 nm blue for acne-related bacteria — while a near-infrared layer (830 nm + 1072 nm, including 6 dedicated 1072 nm collagen diodes) carries energy into the deeper dermis, where visible light stops. Food-grade silicone keeps every LED at the validated 1 cm irradiation distance — the gold standard for LED efficacy. One session: ten minutes, hands-free. Visible change in tone and firmness within 4–8 weeks of regular use.



- ↑ 35% — collagen synthesis
- ↓ 60% — inflammatory redness
- ↑ 28% — skin firmness & elasticity
- ↓ 70% — c. acnes viability



NEAR INFRARED
LIGHT THERAPY



RED LIGHT



BLUE LIGHT



YELLOW LIGHT



PINK LIGHT

SPECTRA, BEEA SKIN, LASER MASK

SPECTRA pairs LED photobiomodulation with true laser precision. 80 quad-chip SMD LEDs cover the full face across four clinically mapped wavelengths — 630 nm red for collagen, 605 nm amber for tone and microcirculation, 460 nm blue for acne-related bacteria and 850 nm near-infrared for deep-tissue circulation — while 38 laser diodes at 660 nm deliver coherent light at the optimal LLLT wavelength, beyond the scattering limits of standard LED masks. Professional-grade irradiance of 35–50 mW/cm². Wireless, with a 2000 mAh battery. Clinical-grade light therapy — worn at home.



collagen activation

acne and sebum control

firmness and fine lines

circulation and recovery

even tone and redness

full-face precision coverage



LOW LEVEL LASER
TECHNOLOGY




RED LIGHT



BLUE LIGHT



AMBER LIGHT



NEAR INFRARED
LIGHT THERAPY

FENIKS, BEEA SKIN, HAIR REGROWTH HELMET




FENIKS treats the scalp as a biological system — not cosmetic territory. 54 VCSEL lasers at 660 nm, the optimal LLLT wavelength for follicle-level photobiomodulation, combine with 70 red LED diodes into a 124-source light field. 88 transparent light-guiding teeth part the hair and carry energy directly to the follicle matrix — eliminating the energy loss standard helmets suffer through hair. Light reaches where growth begins: the follicle extends its growth (anagen) phase, and density becomes visible with consistent use. At 345 g with a flexible modular structure and built-in battery, FENIKS is worn — not performed.

+30% effective irradiation coverage vs standard helmet design

124 light sources — 54 VCSEL + 70 LED combined field

660 nm — optimal LLLT wavelength

ultralight, 40% lighter than rigid helmets



FOLLICLE
ACTIVATION



RED LIGHT




RED LIGHT



FACE GYM II, BEEA SKIN, LIFT & CONTOUR

Facial contours are held by more than 40 muscles — and muscles respond to training. FACE GYM II delivers targeted EMS impulses that contract facial muscle fibres, building tone and structural support beneath the skin. Controlled heat prepares the tissue and boosts circulation for more effective activation, while three LED wavelengths work on the skin itself: red for collagen activity, green for balanced complexion, blue for acne-related bacteria. Three targeted programs, four intensity levels each — and an ergonomic extended handle for precise control along the jawline, cheeks, eyes and lips. Improved facial muscle tone typically within 3–4 weeks of consistent use.


- lifts and defines facial contours
- smooths the look of fine lines
- improves muscle tone and firmness
- supports firmer-looking skin
- enhances jawline definition
- boosts microcirculation
- refines skin texture
- reduces the look of sagging



EMS TECHNOLOGY



HEATING FUNCTION



DUAL VIBRATION
MODES



LED LIGHT



ACTIVE MUSCLE STIMULATION

Delivers targeted electrical impulses that contract facial muscles




FACIAL MUSCLE TONING PROTOCOL

*Stimulates muscle fibers
Supports improved muscle density
Enhances structural support beneath the skin*

WHAT USERS CAN EXPECT WITH REGULAR USE*

- Improved facial muscle tone within 3–4 weeks
- More defined jawline and cheek contours
- Visibly firmer skin appearance
- Reduced look of sagging in lower face
- Enhanced structural lift effect

**Results depend on consistency and individual muscle response.*



THERMAL PERFORMANCE BOOST

Controlled heat enhances circulation and prepares tissue for more effective muscle activation



RED LED COLLAGEN SUPPORT

Promotes collagen activity, helping smooth fine lines and improve skin texture

GREEN LED COMPLEXION BALANCE

Green light supports more even-looking tone and contributes to refined skin radiance





BEEA SKIN

SKIN THAT PERFORMS.

Science-led skincare. Engineered for results — with or without a device.

Technology brought you here. Science keeps you here. BEEA SKIN skincare — formulated for results that last beyond the treatment.

1. COLLAGEN HYDRO CONDUCTIVE BARRIER MASK PRO

Your skin's reset button.

Ceramides rebuild. Panthenol calms. Your skin glows — whether it just had a treatment or just had a Tuesday. Device-compatible · Barrier-restoring · Instant luminosity

2. CONDUCTIVE LUXE GLIDE

The glide that changed the standard.

Zero drag. Zero residue. Zero guessing. Precision conductivity for RF, EMS and microcurrent — with a finish your skin will actually thank you for. Professional grade · Skin-respecting · Clinical control

3. AURA DEEP INFUSION ESSENCE

Peptides that go where serums can't.

Optimised for electroporation. Built for skin that demands more. Matrix renewal, deep hydration and visible calm — in one milky infusion. Electroporation-ready · Peptide complex · Sensitive skin approved

4. VITA SCALP DENSITY SERUM

Stronger follicles. Visible density. Daily.

Lightweight. Fast-absorbing. LLLT-compatible. Because hair longevity starts at the scalp — not the strand. Follicle support · Scalp-balancing · Non-greasy daily serum

Formulated for technology. Trusted by skin.



Marketing



INSTAGRAM

Over the past year, we’ve tested various types of content and video formats to discover what resonates best with our audience. We tried it all – from dynamic Reels and behind-the-scenes clips to interactive polls and informative posts. After analyzing what most effectively increased reach and engagement, we now have a clear direction for shaping our future strategy.

This year on Instagram, our main focus will be on high-quality collaborations with creators and brands that share our values, as well as more informative content that delivers useful and engaging insights. We’re also prioritizing active interaction through conversations, interviews, and Q&As.

We aim to involve our community even more through polls and surveys, fun “tag a friend” posts and videos, and giveaways to show appreciation for our followers. In addition, we’ll be using sponsored content strategically to reach a broader audience and continue building a strong community around our brand.

This is a year of growth, genuine connection, and even better content. Follow us, join the conversation, and be a part of our story!



TIKTOK

Over the past year, we tested various types of content on TikTok to discover what works best and brings the highest reach and engagement. We explored different formats, analyzed audience reactions, and now we know which approach delivers the strongest results.

This year, we’re focusing on dynamic, fun, and easily digestible content. Fast-paced videos that grab attention instantly, strong hooks that keep viewers watching from the very first second, and funny, relatable moments that encourage interaction will be at the core of our TikTok strategy.

We want to create content that entertains, connects, and becomes part of viral trends—while staying authentic and true to our story.

Marketing Timeline 2026

Q2 FAZA - FOUNDATION, PROOF & HYPE START

ROLE OF THE QUARTER

Build the foundation for the year through the SCULPTA hero campaign and begin the cosmetic narrative without a reveal.

Quarter focus

- Rems and Revive as hero devices and proof driver
- market education (RF, collagen, lifting mechanism)
- building trust through a 30-day protocol
- beginning of the cosmetic narrative (device + formula synergy)
- forming an EU influencer network

Hero products

- REMS
- Revive Eye Wand
- conductive media (as an introduction to cosmetics)

Events & collaborations

- BEEA Skin Experience (mini events and workshops)
- collaborations with:
 - aesthetic clinics
 - beauty influencers
- PR & editorials (beauty + wellness)
- Estetica.hr awards 2026

Retail & partner role

- educational webshop (clinical logic)
- landing pages (HR + ENG)
- email onboarding:
- protocols
- education
- laying the foundation for EU distribution

Cosmetic thread (without reveal)

- *education:*
 - skin barrier
 - conductivity
 - permeation
- *questions that open the problem:*
 - “What if your serum doesn’t work with the device?”
- *format:*
 - IG carousel (science)
 - TikTok education
 - blog + SEO articles
- *visuals:* laboratories, textures, microphotography (without product)

Influenceri & content

- European influencers
- focus:
- skincare-educated audience
 - educational content
- format:
- proof content
 - routines
 - education

Campaigns & visibility

- “Use. Define. See the difference.”
 - mechanism-first communication (how the therapies work)
 - before/after as core asset
- educational content:
- why RF for lifting
 - how it affects collagen
 - demo videos (real usage + skin texture)

Marketing Timeline 2026

Q3 FAZA - LAUNCH SEQUENCE & SYSTEM BUILD

ROLE OF THE QUARTER

Build the perception of BEEA as a complete system (device + cosmetics) through three key moments.

Quarter focus

- sequential launch: Face Gym II and Cosmetics (Luxe Glide)
- IFA Berlin (B2B credibility)
- intensive education + proof
- building hype and anticipation

Hero products

- Face Gym II (July launch)
- Conductive Luxe Glide (September launch)

Events & Fairs

- BEEA workshops (device education + protocols)
- Sunny brunch event (lifestyle + skincare experience)
- BEEA summer party (community + brand exposure)
- influencer gathering (Face Gym II + cosmetic tease)
- IFA Berlin (B2B + EU visibility)

Campaigns & visibility

1. JULY – Face Gym II

- Your at-home facial workout.”
- EMS + microcurrent education
- 30-day proof content
- cosmetic hint in all content

2. AUGUST – Cosmetic tease

- visuals:
 - textures
 - packaging without text
- countdown:
 - 0. → 7. → 14. → 30. day
- influencers:
 - press kit, embargo, reaction content

3. SEPTEMBER – Luxe Glide launch

- “The gel engineered for your device.”
- full reveal:
 - formula
 - conductivity
 - performance
- hero video (device + gel difference)
- bundles with devices
- early access (waitlist)
- IFA Berlin as a B2B showcase

Influencers & content

- coordinated launch
- UGC + expert content

format:

- demo
- before/after
- education

simultaneous posts on launch day

Retail & partner role

- bundle strategy (device + gel)
- protocol communication (not single product)
- EU distribution (initial entry)
- remarketing + upsell

Marketing Timeline 2026

Q4 FAZA - CLINICAL LUXURY & GIFTING POWER

ROLE OF THE QUARTER

Position BEEA as a premium alternative to clinical treatments and maximize sales through the gifting season.

Quarter focus

- Laser Mask as a flagship launch
- “clinic at home” positioning
- gifting as the key sales moment
- bundles as the core strategy
- premium perception without aggressive discounts

Hero products

- Laser Mask (main launch)
- Luxe Glide (protocol product)
- Barrier Mask (conditional)
- bundle sets (complete protocol)

Events & collaborations

- BEEA Skin Xmas event (media, influencers, VIP customers)
- BEEA workshops (education + device protocols)
- PR presence (gift guides, beauty editorials, expert articles)
- selective influencer & luxury event collaborations
- retail activations (live demo, experience moments)

Campaigns & visibility

Laser Maska launch

- “Clinical-grade light therapy. Now yours.”
- education:
 - LLLT vs standard LED
 - clinical treatment vs device
- cost-per-use communication
- dermatological / aesthetic angle

Gifting campaigns

- “Clinic at home.”
- “The treatment you can gift.”
- focus:
 - self-care
 - premium gift
- unboxing experience as content

Black Friday

- without aggressive discounts
- bundle > discount
- early access / waitlist
- retargeting burst

Influencers & content

- transformation storytelling
- “my skin journey”
- clinic vs at-home comparisons
- gifting & lifestyle content

Retail & partner role

- premium display (luxury positioning)
- gift-ready packaging
- bundle logic in retail
- focus on the experience, not the price

Fairs

- Beautyworld Dubai (6.–8.10.)
- luxury market
- B2B distribution
- gifting positioning



BEEA SKIN

Our Ambassadors

We believe that in today's world, investing in top influencers and well-known personalities is essential to staying in the spotlight and always ahead of the curve.




@dizdar_el
95.4k followers



@ivana.pirizovich
25.7k followers



@merigoldasic
128k followers




@ania_ddd
137k followers




@korana_marovic
17,1k followers




@lanaklingor
157k followers



@tara_osmic
13.1k followers




@anjaozanic
2,5k followers




@nolaculjak
2,3k followers




@belamakeuplada
61,1k followers




@lunavalas
12,3k followers



@isabellae
461k followers



@miamatic
17,7k followers



@matejasucevic
10k followers





BEEA SKIN

[O NAMA](#) [SHOP](#) [INSPIRACIJA](#) [BEEA STORE](#) [ERIDAN](#) [KONTAKT](#)



Budućnost vaše kože,
u vašim je rukama.

Premium uređaji i profesionalni
rezultati u udobnosti vašeg doma.

BEEA® Community



Rems. BEEA SKIN. anti-ageing uređaj



RF TECHNOLOGY

RF TECHNOLOGY

Radiofrequency heats the tissue, stimulates collagen and elastin synthesis, and improves skin tone and elasticity. It helps reduce wrinkles, provides a lifting effect, enhances glow, and long-termly restores the skin structure. For home use, it is best paired with gels, such as ultrasound gel or aloe vera gel for optimal results.



EMS TECHNOLOGY

EMS TECHNOLOGY

EMS therapy stimulates facial muscles using pulsating currents, encouraging their contraction and strengthening. The result is improved muscle tone, a lifting effect, wrinkle reduction, better circulation, and lymphatic drainage, leading to finer and revitalized skin.



LED THERAPY

LED THERAPY

Light therapy involves the application of different light spectra and colors to address specific facial concerns. LED therapy can have an anti-ageing effect, even out skin tone, reduce inflammation, or soothe the skin. Regular use of LED therapy promotes healthier skin and enhances its natural glow.

BEEA SKIN



EMS TECHNOLOGY



LED THERAPY

The REMS beauty device combines three advanced technologies for skin rejuvenation and care. It is easy to use, suitable for all skin types, and delivers visible results quickly. It provides finer and more radiant skin from the comfort of your home.

BEEA[®]SKIN

Hvala!